

A new clever palm-sized pulse oximeter is launched to suit your needs. MD-680P, a small and compact unit with ergonomic design, is easy carried for spot checking and short-term monitoring. It presents an optimal performance with easy operation. Furthermore, the clear LED digital display brings you an outstanding observation. Those unique and important features make MD-680P become the most necessary and important device in clinical, ambulance and homecare application.

Fast: Catch SpO2 & Pulse rate requires only 6 seconds

Smart: Audible sound assists users to detect patient's status

Accurate: Accuracy ±1% between 100% to 80%

Convenient : Perfect compact pocket size easier for measurement













MD 680P SPECIFICATION

SPECIFICATION

SpO₂

SpO₂ Accuracy: 70% ~ 100% ± 2%

50% ~ 69% ± 3%

< 50% unspecified

Pulse Accuracy: $20 \sim 250 \text{ bpm } \pm 1\%$

Alarm

SpO₂ Alarm Range: 0~100% Pulse Alarm Range 20 ~ 300 bpm

Sound: Audible and Mute selectable

(sound varies on pulse rate)

SpO₂ Probe: Built in

Quick Response: 6-8 seconds

(may vary from pulse rate)

0°℃ ~55°℃ **Temperature:**

(operation & storage)

Humidity: 15% ~ 95% (operation)

10% ~ 95% (storage)

Power Requirement: 3 pcs AAA size batteries

Battery Life: Approximate 16 hours of

continuous use mode

Auto Off: Whenever no measurement

detected after 10 seconds

Dimensions: 90(L) x 62(W) x 22(H) mm

90 g (without battery) Weight:

Quality System: ISO 9001, ISO 13485, GMP

FDA (€ 0434

EN-60601-1, EN-60601-1-1 Safety:

EN-60601-1-2


FEATURES

- Palm sized design for easy hand carry
- Quick measurement of a convenient use
- Light weight deisgn with large LED display
- Easily visible LED indicator during night time
- Quick response of measurement
- Built in SpO₂ probe
- SpO₂ and pulse signal strength indicator
- Low battery indicator
- Audible and mute selectable
- Low power consumption with auto off fuction

STANDARD ACCESSORY

- 3 pcs AAA size alkaline batteries
- Carry bag 2.
- 3. Owner's manual

PANEL DISPLAY



COMDEK INDUSTRIAL CORP.

7Fl.,No1 Alley 6, Lane 235,PaoChiao Rd.,

Hsin Tien, Taipei 23145, Taiwan.

Tel: 886-2-29181266 Fax: 886-2-29181277

http://www.comdek.com Email: medical@comdek.com

DICTRIBUTED DV

DISTRIBUTED BY		

^{*} Specifications are subject to change without notice.